

## PRODUCT CATALOGUE

# Shrinkage Solution

## - from Expanding to Hatched Blastocyst -

- Shrinkage of expanding / expanded blastocysts and hatched blastocysts prior to vitrification procedure contributes to effective vitrification.
- HEPES buffered media for use in the ambient atmosphere. Contains Gentamicin as antibiotic substance.



REF	Code	Details	Contents
91366	VT525-4	Shrinkage Solution	1.0 mL×3

### COMPONENTS

HEPES within Basic culture medium / Sucrose /  
Hydroxypropyl cellulose / Gentamicin

### QUALITY CONTROL


pH 7.2-7.6 / Osmolality / Endotoxin <0.25EU/mL /  
Mouse embryo assay ≥80% / Sterility test (bacteria)

Storage : 2-8°C

Shelf life : 12 months

## RESULTS

It is reported that embryos with a zona pellucida inner diameter >140 μm are more likely to have insufficient replacement of water and cryoprotectant and that artificial shrinkage before vitrification improves implantation, pregnancy, and birth rates.



Reference: J Assist Reprod Genet (2016) 33:461-466

Specification may change without pre-notice for purpose of product improvement.


## Kitazato Corporation

HEADQUATER : 100-10 Yanagishima, Fuji, Shizuoka 416-0932 JAPAN

TOKYO : 1-1-8 Shibadaimon, Minato-ku, Tokyo 105-0012 JAPAN

[Mail contact@kitazato.co.jp](mailto:contact@kitazato.co.jp)

## Applicability Guide for Shrinkage Solution




Artificial shrinkage with a Shrinkage Solution (VT525-4) is recommended.

From expanding blastocysts

## PROTOCOL

- 01** Bring Shrinkage Solution to room temperature and dispense 300 $\mu$ L on Repro Plate.  
The procedure is performed at room temperature.
- 02** Transfer embryo to shrinkage solution and observe for 30 seconds.  
If shrinkage is observed within 30 seconds, proceed to vitrification procedure, specifically ES equilibration.  
If shrinkage is not observed, extend the immerse time or pipette to encourage shrinkage. Proceed to vitrification procedure, specifically ES equilibration when shrinkage is confirmed.

Shrinkage Solution 300 $\mu$ L



## REFERENCES

- J Assist Reprod Genet (2016) 33:461–466