

**SUPPORTED AND
FULLY ECHOGENIC
ET CATHETER**

RYDEN™

Experience the difference.
Clarity and control every time



- **Fully echogenic** catheter with echo-tip on the guide
- **Outer metallic support** for more control
- **Bulb tip** for atraumatic transfer



Provided by
Atsumi Yoshida, MD, PhD,
Kiba Park Clinic, JAPAN

Type	Order Number	Description	Diameter (Fr)	Length (cm)
RYDEN™ Catheter	213312	Supported Catheter with Guide	3	22
RYDEN™ Catheter	213355	Supported Catheter with Guide	3	25
RYDEN™ Catheter	223312	Supported Catheter and Guide with Obturator	3	22
RYDEN™ Catheter	223355	Supported Catheter and Guide with Obturator	3	25
Stylet	244017	Stylet for guide with 22cm catheter	-	-
Stylet	244020	Stylet for guide with 25cm catheter	-	-



Catheter with metallic support that generates rigidity for a better control of the catheter

Echogenicity along the whole catheter, which will provide better visualization along the way to the uterine cavity



Bulb-tip with Echogenic ring that enhances safe, smooth and atraumatic pass-way through the cervix

Pre-curved guide with memory to give smooth pass along the patient's anatomy



Large and user-friendly opening in the handle for insertion of the preloaded catheter into the guide