

Cell Culture

Gamete Buffer

- For oocyte retrieval and flushing.
- For handling and manipulation of oocytes/embryos.
- Suitable to use in ambient atmosphere.



Order Number	Code	Details	Content
93465	GBS-100	Gamete Buffer with Human Serum Albumin	100mL
93466	GBS-500		500mL

COMPONENTS

Calcium Chloride / Essential amino acids / Gentamicin / Glucose / HEPES / Magnesium Sulfate / Non-essential amino acids / Potassium Chloride / Potassium Phosphate / Sodium Bicarbonate / Sodium Chloride / Sodium Lactate / Sodium Pyruvate / Taurine / Human Serum Albumin

QUALITY CONTROL

pH 7.2-7.6 / Osmolality 270-295 mOsm/kg / Endotoxin<0.25EU/mL / Mouse Embryo Assay ≥80% / Sterile Filtration / Sterility Test

Storage: 2-8°C / Shelf Life: 8 months