

Cell Culture

# Hypure® Oil

- Light and Heavy, 2 types are available.
- High-purity paraffin oil without RNA detection.
- Prepared oil which is washed twice with ultrapure water.
- Light shielding glass is used to prevent deterioration of quality by light. Air tightness is kept in the bottle container.



	REF	Code	Contents
Hypure® Oil Heavy	93621	MOHV-100	100mL
	93620	MOHV-500	500mL
Hypure® Oil Light	93521	MOLH-100	100mL
	93520	MOLH-500	500mL

## COMPONENTS

High-purity Paraffin Oil

Specific Gravity    **Light**    **Heavy**  
 0.838g/mL, 25°C    0.862g/mL, 25°C

## QUALITY CONTROL

POV ≤ 0.1 meq/kg / Readily Carbonizable Substances /  
 Endotoxin < 0.25 EU/mL / Mouse Embryo Assay ≥ 80% / Sterility test  
 Storage: 8-30°C / Shelf Life: 18 months