

Micro Manipulation

PBS

Phosphate buffered saline pH 7.4

- Without calcium and magnesium.
- Suitable for washing cells, especially for washing oocytes.
- Suitable to use under ambient atmosphere.



REF	Code	Contents
93352	PBSF-500	500mL
93353	PBSF-100	100mL
93354	PBSFH-500	500mL
93355	PBSFH-100	100mL

*PBSFH is with heparin.

COMPONENTS

Di-sodium Hydrogen Phosphate Anhydrous / Potassium Chloride / Potassium Phosphate / Sodium Chloride
 PBSFH: Heparin

QUALITY CONTROL

pH 7.2-7.6 / Osmolarity 270-295mOsm/L / Endotoxin <0.25EU/mL / Mouse Embryo Assay ≥80% / Sterility Test

Storage: 15-30°C / Shelf Life: 2 years



www.kitazato.co.jp

Manufactured by:
Kitazato Corporation
 100-10 Yanagishima, Fuji, Shizuoka 416-0932 JAPAN

Mail contact@kitazato.co.jp

- Specification may change without pre-notice for purpose of product improvement.
- Kitazato Logo is trade mark of KITAZATO.