


Micro Tools



MicroTools for ICSI

High quality pipettes enable accurate and precise manipulation for ART. Pipettes co-developed and designed by major clinics in Japan ensure high performance in ease of manipulation and fertility rate.



● ICSI Injection Pipette



It is easy to control the sperm with long arm. Short spike minimize damage to cell cytoplasm. Appearances of respective pipettes are inspected one by one with image system.



| | Order Number | Code | ID (μm) | OD (μm) | Contents |
|------------------|--------------|----------|---------|---------|-------------|
| Injection | 71000 | MT-INJ30 | 3.5 | 6 | 10pcs / box |

● ICSI Holding Pipette

Manufactured with special process in order to hold the oocytes flat and horizontal.



| | Order Number | Code | ID (μm) | OD (μm) | Contents |
|----------------|--------------|---------|---------|---------|-------------|
| Holding | 71100 | MT-HD30 | 15 | 100 | 10pcs / box |

*The angle is 30° for all pipettes.