


PRODUCT CATALOGUE

Micro Tools for ICSI

- High quality pipettes enable accurate and precise manipulation for ART.
- Pipettes designed with major clinics in Japan ensure high performance in ease of manipulation and high fertility rate.






ICSI Injection Pipette

Its long arm facilitates better control of sperm and its short spike minimizes damage to cell cytoplasm. Appearances of respective pipettes are inspected one by one with image system.

| REF | Code | ID (μm) | OD (μm) | Contents |
|-------|----------|---------|---------|-------------|
| 71000 | MT-INJ30 | 4 | 6 | 10pcs / box |

MT-INJ30






ICSI Holding Pipette

Manufactured with special process in order to hold the oocytes flat and horizontal.

| REF | Code | ID (μm) | OD (μm) | Contents |
|-------|---------|---------|---------|-------------|
| 71100 | MT-HD30 | 15 | 100 | 10pcs / box |

MT-HD30



Specification may change without pre-notice for purpose of product improvement.

Kitazato Corporation

HEADQUARTER : 100-10 Yanagishima, Fuji, Shizuoka 416-0932 JAPAN

TOKYO : 1-1-8 Shibadaimon, Minato-ku, Tokyo 105-0012 JAPAN

[Mail contact@kitazato.co.jp](mailto:contact@kitazato.co.jp)